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- Walker DK, Brendzel AM, Scotten LN. The new St. Jude Medical Regent mechanical heart valve: laboratory measurements of hydrodynamic performance. *J Heart Valve Dis.* 1999 Nov;8(6):687-696.
- Product and Design Evaluation Center, St. Jude Medical, Inc.
- St. Jude Medical, Inc., Pre-Market Approval Application—Summary of Safety and Effectiveness. SJM Regent heart valve, P810002/S57.
- Pibarot P, Dumesnil JG. Hemodynamic and clinical impact of prosthesis-patient mismatch in the aortic valve position and its prevention. *J Am Coll Cardiol.* 2000 Oct;36(4):1131-1141.
- De Paulis R, Sommariva L, Colagrande L, et al. Regression of left ventricular hypertrophy after aortic valve replacement for aortic stenosis with different valve substitutes. *J Thorac Cardiovasc Surg.* 1998 Oct;116(4):590-598.
- Jin XY, Gibson DG, Yacoub MH, Pepper JR. Perioperative assessment of aortic homograft, Toronto stentless valve, and stented valve in the aortic position. *Ann Thorac Surg.* 1995 Aug;60(2 Suppl):S395-401.
- Jin XY, Zhang ZM, Gibson DG, Yacoub MH, Pepper JR. Effects of valve substitute on changes in left ventricular function and hypertrophy after aortic valve replacement. *Ann Thorac Surg.* 1996 Sep;62(3):683-690.
- Del Rizzo DF, Goldman BS, Christakis GT, David TE. Hemodynamic benefits of the Toronto Stentless Valve. *J Thorac Cardiovasc Surg.* 1996 Dec;112(6):1431-1445; discussion 1445-1446.
- Lund O. Valve replacement for aortic stenosis: the curative potential of early operation. *Scand J Thorac Cardiovasc Surg.* 1993;40(Suppl):1-137.
- He GW, Grunkemeier GL, Gately HL, Furnary AP, Starr A. Up to thirty-year survival after aortic valve replacement in the small aortic root. *Ann Thorac Surg.* 1995 May;59(5):1056-1062.
- Kannel WB, Cobb J. Left ventricular hypertrophy and mortality—results from the Framingham Study. *Cardiology.* 1992;81(4-5):291-298.
- Levy D, Garrison RJ, Savage DD, Kannel WB, Castelli WP. Prognostic implications of echocardiographically determined left ventricular mass in the Framingham Heart Study. *N Engl J Med.* 1990 Dec 13;323(24):1706-1707.
- Dunn FG, McLenachan J, Isles CG, et al. Left ventricular hypertrophy and mortality in hypertension: an analysis of data from the Glasgow Blood Pressure Clinic. *J Hypertens.* 1990 Aug;8(8):775-782.
- Krayenbuehl HP, Hess OM, Monrad ES, Schneider J, Mall G, Turina M. Left ventricular myocardial structure in aortic valve disease before, intermediate, and late after aortic valve replacement. *Circulation.* 1989 Apr;79(4):744-755.
- Orsinelli DA, Aurigemma GP, Battista S, Krendel S, Gaasch WH. Left ventricular hypertrophy and mortality after aortic valve replacement for aortic stenosis. A high risk subgroup identified by preoperative relative wall thickness. *J Am Coll Cardiol.* 1993 Nov 15;22(6):1679-1683.
- Villari B, Vassalli G, Monrad ES, Chiariello M, Turina M, Hess OM. Normalization of diastolic dysfunction in aortic stenosis late after valve replacement. *Circulation.* 1995 May 1;91(9):2353-2358.
- Christakis GT, Joyner CD, Morgan CD, et al. Left ventricular mass regression early after aortic valve replacement. *Ann Thorac Surg.* 1996 Oct;62(4):1084-1089.
- Bach DS, Goldberg M, Sakwa MP, Petracek M, Errett L, Mohr F. Hemodynamics and early performance of the St. Jude Medical Regent aortic valve prosthesis. *J Heart Valve Dis.* 2001 Jul;10(4):436-442.

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