




Diese Zeichnung gilt auch für die entsprechende unsterile Variante!  
 This drawing is valid also for the unsterile set!



Tuch 1000x1000



Verpackung:  
 Tüte 335x500

		Date	IC	Date	Name			
				16.12.2005	R. Sobota	 <p>cardiopulmonary solutions</p> <p>This document contains confidential and proprietary information that is protected by copyright. All design, manufacturing, reproduction, disclosure, use and sale rights associated with this document are reserved by MEDOS AG.</p>		
/	Pos. 2 Linie v. 330 cm auf 250 cm geändert. Siehe Freigabe v. 30.08.2011	08.09.2011	AGI					
						description DE	description GB	
Name v. Medante auf Implanta geändert.		02.09.2011	AGI	Blutkardioplegie-Set 4:1 Modell Implanta		Boodcardioplegia - Set 4:1 Model Implanta		
2	Pos. 2 v. PVC 1/4x3/32, 130 auf PVC 1/4x1/16 330 cm geänd.	30.08.2011	AGI			set-/ drawingnumber		
1	Modell von Raymed in Medante geändert.	14.04.2010	RS			MEH32137_2		
							page	1 / 1